



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: A01430H **Lot #:** 1C07818

Description: C Reactive Protein

Source: Human Fluids

Format: Purified, Liquid

Purification: 99.29% (SDS-PAGE)
Purity: Single band at ~21 kD.

Concentration: 2.82 mg/mL (OD280nm, $E^{0.1\%} = 1.7$)

Buffer: 20 mM Tris-HCl, pH 7.9 to 8.1 containing 280 mM Sodium Chloride, 5 mM Calcium Chloride.

Preservative: 0.09% Sodium Azide

Applications: Specific methodologies have not been tested using this product.

Storage: Store at 2-8°C.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.
No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.

Brenda Dunn

19Mar2018

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY