

Important Note:

Safety Note(s):

5171 Wilfong Road Memphis, TN 38134 USA

Telephone: 901-382-8716 Fax: 901-333-8223

Email: info@meridianlifescience.com www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #:	A82605P	Lot #:	10H21913
Description:	Estriol-6 Horseradish Peroxidase Conjugated Estriol linked at the 6 position.		
Source:	Not Applicable		
Format:	HRP, Liquid Note: Use of sodium azide as a preservative will substantially inhibit the peroxidase.	ne enzyme activity of	horseradish
Concentration:	Not Determined		
Buffer:	Tris containing protein stabilizers.		
Preservative:	0.005% Thimerosal and other additives.		
Applications:	Specific methodologies have not been tested using this product.		
Storage:	Store at 2-8°C.		
Inactivation:	Not Applicable		

CorryThomas

23 March 2017

Signature Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Refer to the appropriate Safety Data Sheet (SDS) for additional information.