

SAFETY DATA SHEET

according to 1907/2006/EC, Article 31

Page 1/4

anti-IgG1 (human), mAb (Biotin)

Revision 0

Revision date 2013-04-02

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

| | |
|--------------|---------------------------------|
| Product name | anti-IgG1 (human), mAb (Biotin) |
| Product code | AG-29B-0007B |

1.3. Details of the supplier of the safety data sheet

| | |
|-----------|--|
| Company | Adipogen International |
| Address | Adipogen AG Schützenstrasse 12 Postfach 335 CH-4410 Liestal |
| Web | www.adipogen.com |
| Telephone | +41-61-926-60-40 |
| Fax | +41-61-926-60-49 |
| Email | info@adipogen.com |

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

| | |
|--------------|-----------------------|
| Main hazards | No Significant Hazard |
|--------------|-----------------------|

2.2. Label elements

| | |
|--------------|------------------------|
| Risk phrases | No Significant Hazard. |
|--------------|------------------------|

SECTION 3: Composition/information on ingredients

3.2. Mixtures

67/548/EEC / 1999/45/EC

| Chemical Name | Index No. | CAS No. | EC No. | REACH Registration Number | Conc. (%w/w) | Classification |
|--|--------------|------------|-----------|---------------------------|--------------|-----------------------|
| Sodium chloride | | 7647-14-5 | | | 0.5 - 1% | |
| Sodium azide (Sodium azide (as NaN ₃)) | 011-004-00-7 | 26628-22-8 | 247-852-1 | | 0 - 0.5% | T+; R28 R32 N; R50/53 |

Description

The product is classified as non hazardous.

SECTION 4: First aid measures

4.1. Description of first aid measures

| | |
|--------------|---|
| Inhalation | May cause irritation to mucous membranes. Move the exposed person to fresh air. |
| Eye contact | May cause irritation to eyes. Rinse immediately with plenty of water for 15 minutes holding the eyelids open. Seek medical attention if irritation or symptoms persist. |
| Skin contact | May cause irritation to skin. Wash off immediately with plenty of soap and water. Remove contaminated clothing. Seek medical attention if irritation or symptoms persist. |
| Ingestion | May cause irritation to mucous membranes. DO NOT INDUCE VOMITING. Seek medical |

anti-IgG1 (human), mAb (Biotin)

Revision 0

Revision date 2013-04-02

4.1. Description of first aid measures

attention if irritation or symptoms persist.

SECTION 5: Firefighting measures**5.1. Extinguishing media**

Use extinguishing media appropriate to the surrounding fire conditions.

5.2. Special hazards arising from the substance or mixture

Burning produces irritating, toxic and obnoxious fumes.

5.3. Advice for firefighters

Wear suitable respiratory equipment when necessary.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures**

Ensure adequate ventilation of the working area.

6.2. Environmental precautions

Do not allow product to enter drains. Prevent further spillage if safe.

6.3. Methods and material for containment and cleaning up

Absorb with inert, absorbent material. Sweep up. Transfer to suitable, labelled containers for disposal. Clean spillage area thoroughly with plenty of water.

SECTION 7: Handling and storage**7.1. Precautions for safe handling**

Avoid contact with eyes and skin. Ensure adequate ventilation of the working area. Adopt best Manual Handling considerations when handling, carrying and dispensing.

7.2. Conditions for safe storage, including any incompatibilities

Keep in a cool, dry, well ventilated area. Keep containers tightly closed. Store in correctly labelled containers.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****8.1.1. Exposure Limit Values**Sodium azide (Sodium azide (as NaN₃))

WEL 8-hr limit ppm: -

WEL 8-hr limit mg/m³: 0.1

WEL 15 min limit ppm: -

WEL 15 min limit mg/m³: 0.3**8.2. Exposure controls****8.2.1. Appropriate engineering controls**

Ensure adequate ventilation of the working area.

8.2.2. Individual protection measures

Wear protective clothing.

Eye / face protection

In case of splashing, wear:.. Approved safety goggles.

Skin protection - Handprotection

Chemical resistant gloves (PVC).

SECTION 9: Physical and chemical properties

anti-IgG1 (human), mAb (Biotin)

Revision 0

Revision date 2013-04-02

9.1. Information on basic physical and chemical properties

| | |
|-------------------|------------------|
| State | Liquid |
| Colour | Colourless |
| Odour | Odourless |
| Solubility | Soluble in water |

SECTION 10: Stability and reactivity

10.2. Chemical stability

| | |
|--|---------------------------------|
| | Stable under normal conditions. |
|--|---------------------------------|

SECTION 11: Toxicological information

11.1. Information on toxicological effects

| | |
|-----------------------|------------------------|
| Acute toxicity | No Significant Hazard. |
|-----------------------|------------------------|

11.1.4. Toxicological Information

| | |
|--|-------------------|
| | No data available |
|--|-------------------|

SECTION 12: Ecological information

12.1. Toxicity

| | |
|--|-------------------|
| | No data available |
|--|-------------------|

Further information

| | |
|--|---|
| | No known adverse environmental effects. |
|--|---|

SECTION 13: Disposal considerations

General information

| | |
|--|---|
| | Dispose of in compliance with all local and national regulations. |
|--|---|

SECTION 14: Transport information

ADR/RID

| | |
|--|--|
| | The product is not classified as dangerous for carriage. |
|--|--|

IMDG

| | |
|--|--|
| | The product is not classified as dangerous for carriage. |
|--|--|

IATA

| | |
|--|--|
| | The product is not classified as dangerous for carriage. |
|--|--|

Further information

| | |
|--|--|
| | The product is not classified as dangerous for carriage. |
|--|--|

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

| | |
|--------------------|--|
| Regulations | COMMISSION REGULATION (EU) No 453/2010 of 20 May 2010 amending Regulation (EC) No 1907/2006 of the European Parliament and of the Council on the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing a European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No 793/93 and Commission Regulation (EC) No 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC. |
|--------------------|--|

SECTION 16: Other information

Other information

| | |
|--|---|
| Text of risk phrases in Section 3 | R28 - Very toxic if swallowed. R32 - Contact with acids liberates very toxic gas. R50/53 - Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment. |
|--|---|

anti-IgG1 (human), mAb (Biotin)

Revision 0

Revision date 2013-04-02

Further information

| | |
|--|--|
| | The information supplied in this Safety Data Sheet is designed only as guidance for the safe use, storage and handling of the product. This information is correct to the best of our knowledge and belief at the date of publication however no guarantee is made to its accuracy. This information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any other process. |
|--|--|
