



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

12/18/2018

Kit Components

Product code	Description
--------------	-------------

D4310	HostZERO Microbial DNA Kit
--------------	-----------------------------------

Components:

D4310-1-20	Host Depletion Solution
D4310-2-5	Microbial Selection Buffer
D4310-3-50	Microbial Selection Enzyme
D3001-2-A	Proteinase K
D3001-2-C	Proteinase K Storage Buffer
R1100-50	DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA...
D4300-1-40	ZymoBIOMICS Lysis Solution
D4300-2-100	ZymoBIOMICS DNA Binding Buffer
D4300-3-50	ZymoBIOMICS DNA Wash Buffer 1
D4300-4-60	ZymoBIOMICS DNA Wash Buffer 2
D4302-5-1, D4302-5-3	ZymoBIOMICS DNase/RNase Free Water



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

1 Identification

- **Product identifier**
- **Trade name:** Host Depletion Solution
- **Article number:** D4310-1-20
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The product is not classified, according to the Globally Harmonized System (GHS).

-
- **Label elements**
 - **GHS label elements** Void
 - **Hazard pictograms** Void
 - **Signal word** Void
 - **Hazard statements** Void
 - **Classification system:**
 - **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Host Depletion Solution

(Contd. of page 1)

3 Composition/information on ingredients

- Chemical characterization: Mixtures
Description: Mixture of the substances listed below with nonhazardous additions.
Dangerous components: Void

4 First-aid measures

- Description of first aid measures
General information: No special measures required.
After inhalation: Supply fresh air; consult doctor in case of complaints.
After skin contact: Generally the product does not irritate the skin.
After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
After swallowing: Do not induce vomiting; immediately call for medical help.
Information for doctor:
Most important symptoms and effects, both acute and delayed No further relevant information available.
Indication of any immediate medical attention and special treatment needed
No further relevant information available.

5 Fire-fighting measures

- Extinguishing media
Suitable extinguishing agents: Use fire fighting measures that suit the environment.
Advice for firefighters
Protective equipment: Wear protective clothing.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures Wear protective clothing.
Environmental precautions: Do not allow to enter sewers/ surface or ground water.
Methods and material for containment and cleaning up: Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Reference to other sections
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
Protective Action Criteria for Chemicals

PAC-1:
None of the ingredients is listed.

PAC-2:
None of the ingredients is listed.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Host Depletion Solution

(Contd. of page 2)

· PAC-3:

None of the ingredients is listed.

7 Handling and storage

- Handling:
· Precautions for safe handling No special measures required.
· Information about protection against explosions and fires: No special measures required.
· Conditions for safe storage, including any incompatibilities
· Storage:
· Requirements to be met by storerooms and receptacles: No special requirements.
· Information about storage in one common storage facility: Not required.
· Further information about storage conditions: None.
· Specific end use(s) Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
· Control parameters
· Components with limit values that require monitoring at the workplace:
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
· Additional information: The lists that were valid during the creation were used as basis.
· Exposure controls
· Personal protective equipment:
· General protective and hygienic measures:
The usual precautionary measures for handling chemicals should be followed.
· Breathing equipment: Not required.
· Protection of hands:
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
· Material of gloves
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
· Penetration time of glove material
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
· Eye protection: Goggles recommended during refilling.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Host Depletion Solution

(Contd. of page 3)

9 Physical and chemical properties

· Information on basic physical and chemical properties

· General Information

· Appearance:

Form:	Liquid
Color:	Light brown
· Odor:	No information available
· Odor threshold:	Not determined.

· pH-value: Not determined.

· Change in condition

Melting point/Melting range:	Undetermined.
Boiling point/Boiling range:	Undetermined.

· Flash point: Not applicable.

· Flammability (solid, gaseous): Not applicable.

· Decomposition temperature: Not determined.

· Auto igniting: Product is not selfigniting.

· Danger of explosion: Product does not present an explosion hazard.

· Explosion limits:

Lower:	Not Applicable
Upper:	Not Applicable

· Vapor pressure: Not determined.

· Density at 20 °C (68 °F): 0.5 g/cm³ (4.1725 lbs/gal)

· Relative density Not determined.

· Vapor density Not determined.

· Evaporation rate Not determined.

· Solubility in / Miscibility with

Water: Not miscible or difficult to mix.

· Partition coefficient (n-octanol/water): Not determined.

· Viscosity:

Dynamic:	Not determined.
Kinematic:	Not determined.

· Solvent content:

VOC content: 0.00 %
0.0 g/l / 0.00 lb/gal

Solids content: 0.0 %

· Other information No further relevant information available.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Host Depletion Solution

(Contd. of page 4)

10 Stability and reactivity

- Reactivity No further relevant information available.
Chemical stability
Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
Possibility of hazardous reactions No dangerous reactions known.
Conditions to avoid No further relevant information available.
Incompatible materials: No further relevant information available.
Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:
Primary irritant effect:
on the skin: No irritant effect.
on the eye: No irritating effect.
Sensitization: No sensitizing effects known.
Additional toxicological information:
The product is not subject to classification according to internally approved calculation methods for preparations:
When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.
Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes:
Water hazard class 2 (Self-assessment): hazardous for water
Do not allow product to reach ground water, water course or sewage system.
Danger to drinking water if even small quantities leak into the ground.

(Contd. on page 6)



Safety Data Sheet acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Host Depletion Solution

(Contd. of page 5)

- Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation: Smaller quantities can be disposed of with household waste.
Uncleaned packagings:
Recommendation: Dispose of container in accordance with local/regional/national and international recommendations.

14 Transport information

Table with 2 columns: Hazard/Requirement and Status. Rows include UN-Number, UN proper shipping name, Transport hazard class(es), Packing group, Environmental hazards, Special precautions for user, Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code, and UN Model Regulation.

15 Regulatory information

- Safety, health and environmental regulations/legislation specific for the substance or mixture
Sara

Section 355 (extremely hazardous substances): None of the ingredients is listed.

Section 313 (Specific toxic chemical listings): None of the ingredients is listed.

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Host Depletion Solution

(Contd. of page 6)

TSCA (Toxic Substances Control Act):

None of the ingredients is listed.

Proposition 65

Chemicals known to cause cancer:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

GHS label elements Void

Hazard pictograms Void

Signal word Void

Hazard statements Void

Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Department issuing SDS:

Zymo Research Corp.
Safety Department
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682

Contact: sds@zymoresearch.com

Date of preparation / last revision 12/18/2018 / -

Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG: International Maritime Code for Dangerous Goods
DOT: US Department of Transportation
IATA: International Air Transport Association

(Contd. on page 8)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 8/8

Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Host Depletion Solution

(Contd. of page 7)

ACGIH: American Conference of Governmental Industrial Hygienists
EINECS: European Inventory of Existing Commercial Chemical Substances
ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
NFPA: National Fire Protection Association (USA)
HMIS: Hazardous Materials Identification System (USA)
VOC: Volatile Organic Compounds (USA, EU)
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
NIOSH: National Institute for Occupational Safety
OSHA: Occupational Safety & Health
TLV: Threshold Limit Value
PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

1 Identification

- **Product identifier**
- **Trade name: Microbial Selection Buffer**
- **Article number:** D4310-2-5
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The product is not classified, according to the Globally Harmonized System (GHS).

-
- **Label elements**
 - **GHS label elements** Void
 - **Hazard pictograms** Void
 - **Signal word** Void
 - **Hazard statements** Void
 - **Classification system:**
 - **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Buffer

(Contd. of page 1)

3 Composition/information on ingredients

- Chemical characterization: Mixtures
Description: Mixture of the substances listed below with nonhazardous additions.
Dangerous components: Void

4 First-aid measures

- Description of first aid measures
General information: No special measures required.
After inhalation: Supply fresh air; consult doctor in case of complaints.
After skin contact: Generally the product does not irritate the skin.
After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
After swallowing: Do not induce vomiting; immediately call for medical help.
Information for doctor:
Most important symptoms and effects, both acute and delayed No further relevant information available.
Indication of any immediate medical attention and special treatment needed
No further relevant information available.

5 Fire-fighting measures

- Extinguishing media
Suitable extinguishing agents: Use fire fighting measures that suit the environment.
Advice for firefighters
Protective equipment: Wear protective clothing.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures Wear protective clothing.
Environmental precautions: Dilute with plenty of water.
Methods and material for containment and cleaning up: Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Reference to other sections
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
Protective Action Criteria for Chemicals

PAC-1:
None of the ingredients is listed.

PAC-2:
None of the ingredients is listed.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Buffer

(Contd. of page 2)

· PAC-3:

None of the ingredients is listed.

7 Handling and storage

- **Handling:**
- **Precautions for safe handling** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s)** Laboratory reagent

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**
The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:** Not required.
- **Protection of hands:**
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Goggles recommended during refilling.

US

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Buffer

(Contd. of page 3)

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information

Appearance:

- Form: Liquid
Color: Clear
Odor: No information available
Odor threshold: Not determined.

pH-value: Not determined.

Change in condition

- Melting point/Melting range: Undetermined.
Boiling point/Boiling range: Undetermined.

Flash point: Not applicable.

Flammability (solid, gaseous): Not applicable.

Decomposition temperature: Not determined.

Auto igniting: Product is not selfigniting.

Danger of explosion: Product does not present an explosion hazard.

Explosion limits:

- Lower: Not Applicable
Upper: Not Applicable

Vapor pressure: Not determined.

Density: Not determined.

Relative density: Not determined.

Vapor density: Not determined.

Evaporation rate: Not determined.

Solubility in / Miscibility with

Water: Fully miscible.

Partition coefficient (n-octanol/water): Not determined.

Viscosity:

- Dynamic: Not determined.
Kinematic: Not determined.

Solvent content:

VOC content: 0.00 %
0.0 g/l / 0.00 lb/gal

Solids content: 0.0 %

Other information: No further relevant information available.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Buffer

(Contd. of page 4)

10 Stability and reactivity

- Reactivity No further relevant information available.
Chemical stability
Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
Possibility of hazardous reactions No dangerous reactions known.
Conditions to avoid No further relevant information available.
Incompatible materials: No further relevant information available.
Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:
Primary irritant effect:
on the skin: No irritant effect.
on the eye: No irritating effect.
Sensitization: No sensitizing effects known.
Additional toxicological information:
The product is not subject to classification according to internally approved calculation methods for preparations:
When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.
Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes: Not known to be hazardous to water.
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.

(Contd. on page 6)



Safety Data Sheet acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Buffer

Other adverse effects No further relevant information available.

(Contd. of page 5)

13 Disposal considerations

- Waste treatment methods
Recommendation: Smaller quantities can be disposed of with household waste.
Uncleaned packagings:
Recommendation: Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Regulatory information and Status. Rows include UN-Number, UN proper shipping name, Transport hazard class(es), Packing group, Environmental hazards, Special precautions for user, Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code, and UN Model Regulation.

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture
Sara

Table with 2 columns: Regulatory section and Status. Rows include Section 355 (extremely hazardous substances), Section 313 (Specific toxic chemical listings), and TSCA (Toxic Substances Control Act).

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Buffer

(Contd. of page 6)

· **Proposition 65**

· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements** Void

· **Hazard pictograms** Void

· **Signal word** Void

· **Hazard statements** Void

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 12/18/2018 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

(Contd. on page 8)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 8/8

Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Buffer

(Contd. of page 7)

NFPA: National Fire Protection Association (USA)
HMIS: Hazardous Materials Identification System (USA)
VOC: Volatile Organic Compounds (USA, EU)
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
NIOSH: National Institute for Occupational Safety
OSHA: Occupational Safety & Health
TLV: Threshold Limit Value
PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018


1 Identification

- **Product identifier**
- **Trade name:** *Microbial Selection Enzyme*
- **Article number:** D4310-3-50
- **CAS Number:**
56-81-5
- **EC number:**
200-289-5
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The substance is not classified, according to the Globally Harmonized System (GHS).

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



Health = 0
Fire = 0
Reactivity = 0

- **HMIS-ratings (scale 0 - 4)**



Health = 0
Fire = 0
Reactivity = 0

- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.

(Contd. on page 2)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Enzyme

· **vPvB:** Not applicable.

(Contd. of page 1)

3 Composition/information on ingredients

- **Chemical characterization: Substances**
- **CAS No. Description**
56-81-5 glycerol
- **Identification number(s)**
- **EC number:** 200-289-5

4 First-aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.
- **After swallowing:** Do not induce vomiting; immediately call for medical help.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Wear protective clothing.
- **Environmental precautions:**
Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Enzyme

(Contd. of page 2)

Protective Action Criteria for Chemicals

Table with 2 columns: Criteria (PAC-1, PAC-2, PAC-3) and Value (45 mg/m³, 180 mg/m³, 1,100 mg/m³)

7 Handling and storage

- Handling:
Precautions for safe handling: No special measures required.
Information about protection against explosions and fires: No special measures required.
Conditions for safe storage, including any incompatibilities
Storage:
Requirements to be met by storerooms and receptacles: No special requirements.
Information about storage in one common storage facility: Not required.
Further information about storage conditions: None.
Specific end use(s): Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
Control parameters

Table with 2 columns: Component (CAS, PEL, TLV) and Value (56-81-5 glycerol, Long-term value: 15* 5** mg/m³ mist; *total dust **respirable fraction, TLV withdrawn-insufficient data human occup. exp.)

- Additional information: The lists that were valid during the creation were used as basis.
Exposure controls
Personal protective equipment:
General protective and hygienic measures: The usual precautionary measures for handling chemicals should be followed.
Breathing equipment: Not required.
Protection of hands: The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Enzyme

(Contd. of page 3)

- Material of gloves
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.
Penetration time of glove material
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
Eye protection: Goggles recommended during refilling.

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information

Appearance:

- Form: Liquid
Color: Clear
Odor: No information available
Odor threshold: Not determined.

pH-value: Not determined.

Change in condition

- Melting point/Melting range: 18.2 °C (64.8 °F)
Boiling point/Boiling range: 290 °C (554 °F)

Flash point: 160 °C (320 °F)

Flammability (solid, gaseous): Not applicable.

Ignition temperature: 400 °C (752 °F)

Decomposition temperature: Not determined.

Auto igniting: Not determined.

Danger of explosion: Product does not present an explosion hazard.

Explosion limits:

- Lower: Not Applicable
Upper: Not Applicable

Vapor pressure at 20 °C (68 °F): <0.1 hPa (<0.1 mm Hg)

Density at 20 °C (68 °F): 1.26 g/cm³ (10.5147 lbs/gal)

Relative density: Not determined.

Vapor density: Not determined.

Evaporation rate: Not determined.

Solubility in / Miscibility with

Water: Fully miscible.

Partition coefficient (n-octanol/water): Not determined.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Enzyme

(Contd. of page 4)

· Viscosity:	
Dynamic:	Not determined.
Kinematic:	Not determined.
Organic solvents:	100.0 %
VOC content:	0.00 %
	0.0 g/l / 0.00 lb/gal
Solids content:	0.0 %
· Other information	No further relevant information available.

10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**

· **LD/LC50 values that are relevant for classification:**

CAS: 56-81-5 glycerol

Oral	LD50	12,600 mg/kg (rat)
------	------	--------------------

- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**
When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.
The substance is not subject to classification.

- **Carcinogenic categories**

· **IARC (International Agency for Research on Cancer)**

Substance is not listed.

· **NTP (National Toxicology Program)**

Substance is not listed.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Enzyme

(Contd. of page 5)

OSHA-Ca (Occupational Safety & Health Administration)

Substance is not listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.
Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes:
Water hazard class 1 (Assessment by list): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation: Smaller quantities can be disposed of with household waste.
Uncleaned packagings:
Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

- UN-Number
DOT, ADN, IMDG, IATA not regulated
UN proper shipping name
DOT, ADN, IMDG, IATA not regulated
Transport hazard class(es)
DOT, ADN, IMDG, IATA
Class not regulated
Packing group
DOT, IMDG, IATA not regulated

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Enzyme

(Contd. of page 6)

- **Environmental hazards:** Not applicable.
- **Special precautions for user** Not applicable.
- **Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code** Not applicable.
- **UN "Model Regulation":** not regulated

15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Sara**

· **Section 355 (extremely hazardous substances):**

Substance is not listed.

· **Section 313 (Specific toxic chemical listings):**

Substance is not listed.

· **TSCA (Toxic Substances Control Act):**

Substance is listed.

· **Proposition 65**

· **Chemicals known to cause cancer:**

Substance is not listed.

· **Chemicals known to cause reproductive toxicity for females:**

Substance is not listed.

· **Chemicals known to cause reproductive toxicity for males:**

Substance is not listed.

· **Chemicals known to cause developmental toxicity:**

Substance is not listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

Substance is not listed.

· **TLV (Threshold Limit Value established by ACGIH)**

Substance is not listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

Substance is not listed.

- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void

(Contd. on page 8)

**Safety Data Sheet**
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/17/2018

Trade name: Microbial Selection Enzyme

(Contd. of page 7)

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.**16 Other information**

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com· **Date of preparation / last revision** 12/18/2018 / -· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/07/2015

1 Identification

- **Product identifier**
- **Trade name: Proteinase K**
- **Article number:** D3001-2-A, D3001-2-B, D3001-2-F
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The product is not classified, according to the Globally Harmonized System (GHS).

-
- **Label elements**
 - **GHS label elements** Void
 - **Hazard pictograms** Void
 - **Signal word** Void
 - **Hazard statements** Void
 - **Classification system:**
 - **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/07/2015

Trade name: Proteinase K

(Contd. of page 1)

3 Composition/information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** Mixture of the substances listed below with nonhazardous additions.
- **Dangerous components:** Void

4 First-aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.
- **After swallowing:** Do not induce vomiting; immediately call for medical help.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Use fire fighting measures that suit the environment.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Wear protective clothing.
- **Environmental precautions:** Dilute with plenty of water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
- **Protective Action Criteria for Chemicals**

· **PAC-1:**

None of the ingredients is listed.

· **PAC-2:**

None of the ingredients is listed.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/07/2015

Trade name: Proteinase K

(Contd. of page 2)

· **PAC-3:**

None of the ingredients is listed.

7 Handling and storage

- **Handling:**
- **Precautions for safe handling** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s)** Laboratory reagent

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**
The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:** Not required.
- **Protection of hands:**
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Goggles recommended during refilling.

US

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/07/2015

Trade name: Proteinase K

(Contd. of page 3)

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information

Appearance:

Form: Crystalline powder

Color: Whitish

Odor: Odorless

Odor threshold: Not determined.

pH-value: Not determined.

Change in condition

Melting point/Melting range: Undetermined.

Boiling point/Boiling range: Undetermined.

Flash point: Not applicable.

Flammability (solid, gaseous): Not applicable.

Decomposition temperature: Not determined.

Auto igniting: Product is not selfigniting.

Danger of explosion: Product does not present an explosion hazard.

Explosion limits:

Lower: Not Applicable

Upper: Not Applicable

Vapor pressure: Not determined.

Density: Not determined.

Relative density: Not determined.

Vapor density: Not determined.

Evaporation rate: Not determined.

Solubility in / Miscibility with

Water: Fully miscible.

Partition coefficient (n-octanol/water): Not determined.

Viscosity:

Dynamic: Not determined.

Kinematic: Not determined.

Solvent content:

VOC content: 0.00 %
0.0 g/l / 0.00 lb/gal

Solids content: 0.0 %

Other information: No further relevant information available.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/07/2015

Trade name: Proteinase K

(Contd. of page 4)

10 Stability and reactivity

- Reactivity No further relevant information available.
Chemical stability
Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
Possibility of hazardous reactions No dangerous reactions known.
Conditions to avoid No further relevant information available.
Incompatible materials: No further relevant information available.
Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:
Primary irritant effect:
on the skin: No irritant effect.
on the eye: No irritating effect.
Sensitization: No sensitizing effects known.
Additional toxicological information:
The product is not subject to classification according to internally approved calculation methods for preparations:
When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.
Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes: Not hazardous for water.
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.

(Contd. on page 6)



Safety Data Sheet acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/07/2015

Trade name: Proteinase K

(Contd. of page 5)

Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation: Smaller quantities can be disposed of with household waste.
Uncleaned packagings:
Recommendation: Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Regulatory requirement and Status. Rows include UN-Number, UN proper shipping name, Transport hazard class(es), Packing group, Environmental hazards, Special precautions for user, Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code, and UN Model Regulation.

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture
Sara

Table with 2 columns: Regulatory section and Status. Rows include Section 355 (extremely hazardous substances), Section 313 (Specific toxic chemical listings), and TSCA (Toxic Substances Control Act).

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/07/2015

Trade name: Proteinase K

(Contd. of page 6)

· **Proposition 65**

· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements** Void

· **Hazard pictograms** Void

· **Signal word** Void

· **Hazard statements** Void

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 12/18/2018 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

(Contd. on page 8)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 8/8

Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 12/07/2015

Trade name: Proteinase K

(Contd. of page 7)

NFPA: National Fire Protection Association (USA)
HMIS: Hazardous Materials Identification System (USA)
VOC: Volatile Organic Compounds (USA, EU)
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
NIOSH: National Institute for Occupational Safety
OSHA: Occupational Safety & Health
TLV: Threshold Limit Value
PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 11/13/2015

1 Identification

- **Product identifier**
- **Trade name: Proteinase K Storage Buffer**
- **Article number:** D3001-2-C, D3001-2-D, D3001-2-E
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The product is not classified, according to the Globally Harmonized System (GHS).

-
- **Label elements**
 - **GHS label elements** Void
 - **Hazard pictograms** Void
 - **Signal word** Void
 - **Hazard statements** Void
 - **Classification system:**
 - **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 11/13/2015

Trade name: Proteinase K Storage Buffer

(Contd. of page 1)

3 Composition/information on ingredients

- Chemical characterization: Mixtures
Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

Table with 2 columns: Component (CAS: 56-81-5 glycerol) and Concentration (50%)

4 First-aid measures

- Description of first aid measures
General information: No special measures required.
After inhalation: Supply fresh air; consult doctor in case of complaints.
After skin contact: Generally the product does not irritate the skin.
After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
After swallowing: Do not induce vomiting; immediately call for medical help.
Information for doctor:
Most important symptoms and effects, both acute and delayed No further relevant information available.
Indication of any immediate medical attention and special treatment needed
No further relevant information available.

5 Fire-fighting measures

- Extinguishing media
Suitable extinguishing agents: Use fire fighting measures that suit the environment.
Advice for firefighters
Protective equipment: Wear protective clothing.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures Wear protective clothing.
Environmental precautions:
Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water.
Methods and material for containment and cleaning up:
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Reference to other sections
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
Protective Action Criteria for Chemicals

PAC-1:

Table with 2 columns: Component (CAS: 56-81-5 glycerol) and Concentration (45 mg/m3)

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 11/13/2015

Trade name: Proteinase K Storage Buffer

(Contd. of page 2)

PAC-2:

Table with 2 columns: CAS: 56-81-5 glycerol, 180 mg/m³

PAC-3:

Table with 2 columns: CAS: 56-81-5 glycerol, 1,100 mg/m³

7 Handling and storage

- Handling:
Precautions for safe handling: No special measures required.
Information about protection against explosions and fires: No special measures required.
Conditions for safe storage, including any incompatibilities
Storage:
Requirements to be met by storerooms and receptacles: No special requirements.
Information about storage in one common storage facility: Not required.
Further information about storage conditions: None.
Specific end use(s): Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
Control parameters

Components with limit values that require monitoring at the workplace:

CAS: 56-81-5 glycerol

PEL Long-term value: 15* 5** mg/m³
mist; *total dust **respirable fraction

TLV TLV withdrawn-insufficient data human occup. exp.

- Additional information: The lists that were valid during the creation were used as basis.
Exposure controls
Personal protective equipment:
General protective and hygienic measures: The usual precautionary measures for handling chemicals should be followed.
Breathing equipment: Not required.
Protection of hands: The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Material of gloves: The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 11/13/2015

Trade name: Proteinase K Storage Buffer

(Contd. of page 3)

- Penetration time of glove material
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
Eye protection: Goggles recommended during refilling.

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information

Appearance:

- Form: Liquid
Color: Clear
Odor: Odorless
Odor threshold: Not determined.

pH-value: Not determined.

Change in condition

- Melting point/Melting range: Undetermined.
Boiling point/Boiling range: 290 °C (554 °F)

Flash point: 160 °C (320 °F)

Flammability (solid, gaseous): Not applicable.

Ignition temperature: 400 °C (752 °F)

Decomposition temperature: Not determined.

Auto igniting: Product is not selfigniting.

Danger of explosion: Product does not present an explosion hazard.

Explosion limits:

- Lower: 0.9 Vol %
Upper: 0.0 Vol %

Vapor pressure at 20 °C (68 °F): 0.1 hPa (0.1 mm Hg)

- Density: Not determined.
Relative density: Not determined.
Vapor density: Not determined.
Evaporation rate: Not determined.

Solubility in / Miscibility with

Water: Fully miscible.

Partition coefficient (n-octanol/water): Not determined.

Viscosity:

- Dynamic: Not determined.
Kinematic: Not determined.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 11/13/2015

Trade name: Proteinase K Storage Buffer

(Contd. of page 4)

Table with solvent content (Organic solvents: 50.0 %, VOC content: 0.00 %, 0.0 g/l / 0.00 lb/gal), solids content (0.0 %), and other information (No further relevant information available).

10 Stability and reactivity

- Reactivity: No further relevant information available.
Chemical stability
Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
Possibility of hazardous reactions: No dangerous reactions known.
Conditions to avoid: No further relevant information available.
Incompatible materials: No further relevant information available.
Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects
Acute toxicity
Primary irritant effect:
on the skin: No irritant effect.
on the eye: No irritating effect.
Sensitization: No sensitizing effects known.
Additional toxicological information:
The product is not subject to classification according to internally approved calculation methods for preparations: When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 11/13/2015

Trade name: Proteinase K Storage Buffer

(Contd. of page 5)

- Persistence and degradability No further relevant information available.
· Behavior in environmental systems:
· Bioaccumulative potential No further relevant information available.
· Mobility in soil No further relevant information available.
· Additional ecological information:
· General notes:
Water hazard class 1 (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
· Results of PBT and vPvB assessment
· PBT: Not applicable.
· vPvB: Not applicable.
· Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
· Recommendation: Smaller quantities can be disposed of with household waste.
· Uncleaned packagings:
· Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.
· Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Hazard/Requirement and Status. Rows include UN-Number, UN proper shipping name, Transport hazard class(es), Packing group, Environmental hazards, Special precautions for user, Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code, and UN Model Regulation.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 11/13/2015

Trade name: Proteinase K Storage Buffer

(Contd. of page 6)

15 Regulatory information

- *Safety, health and environmental regulations/legislation specific for the substance or mixture*
- *Sara*

· **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**

None of the ingredients is listed.

· **TSCA (Toxic Substances Control Act):**

All ingredients are listed.

· **Proposition 65**

· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements** Void

· **Hazard pictograms** Void

· **Signal word** Void

· **Hazard statements** Void

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.
Safety Department
17062 Murphy Ave.
Irvine, CA 92614
USA

(Contd. on page 8)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 8/8

Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 11/13/2015

Trade name: Proteinase K Storage Buffer

(Contd. of page 7)

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 12/18/2018 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 05/06/2016

1 Identification

- **Product identifier**
- **Trade name:** DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Blood Kit, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit
- **Article number:**
R1100-50, R1100-250, R1100-8-S, R1200-25, R1200-125, R1150, R1101, R1102, R1103, R1104, R1105, R1106, R1107, R1108, R1109, R1210
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**



Acute Tox. 4 H302 Harmful if swallowed.
 Skin Irrit. 2 H315 Causes skin irritation.
 Eye Irrit. 2B H320 Causes eye irritation.

- **Label elements**
- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms** GHS07
- **Signal word** Warning
- **Hazard statements**
Harmful if swallowed.
Causes skin and eye irritation.
- **Precautionary statements**
If medical advice is needed, have product container or label at hand.
Keep out of reach of children.
Read label before use.
Wash thoroughly after handling.
Do not eat, drink or smoke when using this product.

(Contd. on page 2)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Blood Kit, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 1)

- Wear protective gloves.
- If swallowed: Call a poison center/doctor if you feel unwell.
- Rinse mouth.
- If on skin: Wash with plenty of water.
- Specific treatment (see on this label).
- If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- If skin irritation occurs: Get medical advice/attention.
- Take off contaminated clothing and wash it before reuse.
- If eye irritation persists: Get medical advice/attention.
- Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



· **HMS-ratings (scale 0 - 4)**



· **Other hazards**

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** This product is a proprietary solution.

· **Dangerous components:** Void

4 First-aid measures

· **Description of first aid measures**

· **After inhalation:** Supply fresh air; consult doctor in case of complaints.

· **After skin contact:** Generally the product does not irritate the skin.

· **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.

· **After swallowing:** Do not induce vomiting; immediately call for medical help.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Blood Kit, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 2)

- Information for doctor:
Most important symptoms and effects, both acute and delayed No further relevant information available.
Indication of any immediate medical attention and special treatment needed No further relevant information available.

5 Fire-fighting measures

- Extinguishing media
Suitable extinguishing agents: Not considered a fire hazard. Use extinguishing media appropriate surrounding fire conditions. Use fire fighting measures that suit the environment.
Advice for firefighters
Protective equipment: Wear protective clothing.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures Wear protective clothing.
Environmental precautions: Dilute with plenty of water.
Methods and material for containment and cleaning up: Ensure adequate ventilation. Dispose contaminated material as waste according to item 13.
Reference to other sections See Section 7 for information on safe handling. See Section 8 for information on personal protection equipment. See Section 13 for disposal information.
Protective Action Criteria for Chemicals

PAC-1: None of the ingredients is listed.

PAC-2: None of the ingredients is listed.

PAC-3: None of the ingredients is listed.

7 Handling and storage

- Handling:
Precautions for safe handling No special precautions are necessary if used correctly.
Information about protection against explosions and fires: No special measures required.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Blood Kit, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 3)

- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep receptacle tightly sealed.
- **Specific end use(s)** Laboratory reagent

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:** Wash hands before breaks and at the end of work.
- **Breathing equipment:** Not required.
- **Protection of hands:** Use appropriate gloves. Avoid contact with skin.
- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Avoid contact with eyes.
- **Body protection:**
Avoid contact with eyes, skin, and clothing. Wash thoroughly after handling. Remove contaminated clothing and wash before use. Avoid ingestion and inhalation.

9 Physical and chemical properties

- **Information on basic physical and chemical properties**
- **General Information**
- **Appearance:**

· Form:	Liquid
· Color:	Colorless
· Odor:	No information available
· Odor threshold:	Not determined.
- **pH-value:** 5.0-7.0

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Blood Kit, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 4)

Table with 2 columns: Property and Value. Rows include: Change in condition (Melting point, Boiling point), Flash point, Flammability, Decomposition temperature, Auto igniting, Danger of explosion, Explosion limits (Lower, Upper), Vapor pressure, Density, Relative density, Vapor density, Evaporation rate, Solubility in / Miscibility with Water, Partition coefficient (n-octanol/water), Viscosity (Dynamic, Kinematic), Solvent content (VOC content), Solids content, and Other information.

10 Stability and reactivity

- Reactivity Strong oxidizing agents and acids.
Possibility of hazardous reactions Strong oxidizing agents and acids.
Conditions to avoid Strong oxidizing agents and acids.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Blood Kit, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

Incompatible materials: Strong oxidizing agents and acids.

(Contd. of page 5)

11 Toxicological information

- Information on toxicological effects
Acute toxicity:
Primary irritant effect:
on the eye: Irritating to eyes.
Additional toxicological information:
Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.
Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects No further relevant information available.

13 Disposal considerations

- Uncleaned packagings:
Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Blood Kit, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 6)

14 Transport information

· <i>UN-Number</i> · <i>DOT, ADN, IMDG, IATA</i>	Not dangerous goods. not regulated
· <i>UN proper shipping name</i> · <i>DOT, ADN, IMDG, IATA</i>	Not dangerous goods. not regulated
· <i>Transport hazard class(es)</i> · <i>DOT, ADN, IMDG, IATA</i> · <i>Class</i>	Not dangerous goods. not regulated
· <i>Packing group</i> · <i>DOT, IMDG, IATA</i>	Not dangerous goods. not regulated
· <i>Environmental hazards:</i>	Not applicable.
· <i>Special precautions for user</i>	Not applicable.
· <i>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</i>	Not applicable.
· <i>UN "Model Regulation":</i>	not regulated

15 Regulatory information

· <i>Safety, health and environmental regulations/legislation specific for the substance or mixture</i> · <i>Sara</i>	
· <i>Section 355 (extremely hazardous substances):</i>	None of the ingredients is listed.
· <i>Section 313 (Specific toxic chemical listings):</i>	None of the ingredients is listed.
· <i>TSCA (Toxic Substances Control Act):</i>	None of the ingredients is listed.
· <i>Proposition 65</i>	
· <i>Chemicals known to cause cancer:</i>	None of the ingredients is listed.
· <i>Chemicals known to cause reproductive toxicity for females:</i>	None of the ingredients is listed.
· <i>Chemicals known to cause reproductive toxicity for males:</i>	None of the ingredients is listed.

(Contd. on page 8)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Blood Kit, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 7)

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

· Carcinogenic categories

· EPA (Environmental Protection Agency)

None of the ingredients is listed.

· TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

· GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

· Hazard pictograms GHS07

· Signal word Warning

· Hazard statements

Harmful if swallowed.

Causes skin and eye irritation.

· Precautionary statements

If medical advice is needed, have product container or label at hand.

Keep out of reach of children.

Read label before use.

Wash thoroughly after handling.

Do not eat, drink or smoke when using this product.

Wear protective gloves.

If swallowed: Call a poison center/doctor if you feel unwell.

Rinse mouth.

If on skin: Wash with plenty of water.

Specific treatment (see on this label).

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

If skin irritation occurs: Get medical advice/attention.

Take off contaminated clothing and wash it before reuse.

If eye irritation persists: Get medical advice/attention.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Department issuing SDS:

Zymo Research Corp.

Safety Department

(Contd. on page 9)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 9/9

Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Blood Kit, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 8)

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 12/18/2018 / -

· **Abbreviations and acronyms:**

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

Acute Tox. 4: Acute toxicity – Category 4

Skin Irrit. 2: Skin corrosion/irritation – Category 2

Eye Irrit. 2B: Serious eye damage/eye irritation – Category 2B

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

1 Identification

- **Product identifier**
- **Trade name:** ZymoBIOMICS Lysis Solution
- **Article number:** D4300-1-40, D4300-1-150
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The product is not classified, according to the Globally Harmonized System (GHS).

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**

HEALTH	0	Health = 0
FIRE	0	Fire = 0
REACTIVITY	0	Reactivity = 0

- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS Lysis Solution

(Contd. of page 1)

3 Composition/information on ingredients

- Chemical characterization: Mixtures
Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

Table with 3 columns: CAS number, Component name, and Concentration. Row 1: CAS: 6381-92-6, Edetate Disodium, Dihydrate, 20%

4 First-aid measures

- Description of first aid measures
General information: No special measures required.
After inhalation: Supply fresh air; consult doctor in case of complaints.
After skin contact: Generally the product does not irritate the skin.
After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
After swallowing: Do not induce vomiting; immediately call for medical help.
Information for doctor:
Most important symptoms and effects, both acute and delayed No further relevant information available.
Indication of any immediate medical attention and special treatment needed No further relevant information available.

5 Fire-fighting measures

- Extinguishing media
Suitable extinguishing agents: CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam. Use fire fighting measures that suit the environment.
Advice for firefighters
Protective equipment: Wear protective clothing.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures Wear protective clothing.
Environmental precautions: Dilute with plenty of water.
Methods and material for containment and cleaning up: Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Reference to other sections: See Section 7 for information on safe handling. See Section 8 for information on personal protection equipment. See Section 13 for disposal information.
Protective Action Criteria for Chemicals

PAC-1:

Table with 3 columns: CAS number, Component name, and Concentration. Row 1: CAS: 6381-92-6, Edetate Disodium, Dihydrate, 30 mg/m3. Row 2: CAS: 77-86-1, trometamol, 18 mg/m3

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS Lysis Solution

(Contd. of page 2)

· PAC-2:		
CAS: 6381-92-6	Edetate Disodium, Dihydrate	330 mg/m ³
CAS: 77-86-1	trometamol	190 mg/m ³
· PAC-3:		
CAS: 6381-92-6	Edetate Disodium, Dihydrate	2,000 mg/m ³
CAS: 77-86-1	trometamol	1,200 mg/m ³

7 Handling and storage

- **Handling:**
- **Precautions for safe handling** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s)** Laboratory reagent

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**
The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:** Not required.
- **Protection of hands:**
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS Lysis Solution

(Contd. of page 3)

- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Goggles recommended during refilling.

9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

- **Form:** Liquid
- **Color:** Clear
- **Odor:** No information available
- **Odor threshold:** Not determined.

· **pH-value:** Not determined.

· **Change in condition**

- **Melting point/Melting range:** Undetermined.
- **Boiling point/Boiling range:** Undetermined.

· **Flash point:** Not applicable.

· **Flammability (solid, gaseous):** Not applicable.

· **Decomposition temperature:** Not determined.

· **Auto igniting:** Product is not selfigniting.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Explosion limits:**

- **Lower:** Not Applicable
- **Upper:** Not Applicable

· **Vapor pressure:** Not determined.

- **Density:** Not determined.
- **Relative density** Not determined.
- **Vapor density** Not determined.
- **Evaporation rate** Not determined.

· **Solubility in / Miscibility with**

· **Water:** Fully miscible.

· **Partition coefficient (n-octanol/water):** Not determined.

· **Viscosity:**

- **Dynamic:** Not determined.
- **Kinematic:** Not determined.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS Lysis Solution

(Contd. of page 4)

Table with 2 columns: Property (Solvent content, VOC content, Solids content, Other information) and Value (0.00 %, 0.0 g/l / 0.00 lb/gal, 0.0 %, No further relevant information available).

10 Stability and reactivity

- Reactivity, Chemical stability, Thermal decomposition / conditions to be avoided, Possibility of hazardous reactions, Conditions to avoid, Incompatible materials, Hazardous decomposition products.

11 Toxicological information

- Information on toxicological effects, Acute toxicity.

LD/LC50 values that are relevant for classification:

CAS: 6381-92-6 Edetate Disodium, Dihydrate

Oral LD50 2,000 mg/kg (rat)

- Primary irritant effect: on the skin, on the eye, Sensitization, Additional toxicological information.

The product is not subject to classification according to internally approved calculation methods for preparations: When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS Lysis Solution

(Contd. of page 5)

12 Ecological information

- Toxicity
· Aquatic toxicity: No further relevant information available.
· Persistence and degradability No further relevant information available.
· Behavior in environmental systems:
· Bioaccumulative potential No further relevant information available.
· Mobility in soil No further relevant information available.
· Additional ecological information:
· General notes: Not hazardous for water.
· Results of PBT and vPvB assessment
· PBT: Not applicable.
· vPvB: Not applicable.
· Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
· Recommendation: Smaller quantities can be disposed of with household waste.
· Uncleaned packagings:
· Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.
· Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Description and Regulation Status. Rows include UN-Number, UN proper shipping name, Transport hazard class(es), Packing group, Environmental hazards, Special precautions for user, and Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code.

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS Lysis Solution

(Contd. of page 6)

UN "Model Regulation": not regulated

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture

Sara

Section 355 (extremely hazardous substances):

None of the ingredients is listed.

Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

TSCA (Toxic Substances Control Act):

CAS: 7647-14-5 Sodium chloride

CAS: 77-86-1 trometamol

TSCA new (21st Century Act) (Substances not listed)

CAS: 6381-92-6 Edetate Disodium, Dihydrate

Proposition 65

Chemicals known to cause cancer:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

GHS label elements Void

Hazard pictograms Void

Signal word Void

Hazard statements Void

Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

US

(Contd. on page 8)

**Safety Data Sheet**
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS Lysis Solution

(Contd. of page 7)

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Department issuing SDS:

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· Contact: sds@zymoresearch.com**· Date of preparation / last revision** 12/18/2018 / -**· Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

1 Identification

- **Product identifier**
- **Trade name:** ZymoBIOMICS DNA Binding Buffer
- **Article number:** D4300-2-100, D4300-2-250
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**



GHS05 Corrosion

Skin Corr. 1C H314 Causes severe skin burns and eye damage.

Eye Dam. 1 H318 Causes serious eye damage.



GHS07

Acute Tox. 4 H302 Harmful if swallowed.

Aquatic Chronic 3 H412 Harmful to aquatic life with long lasting effects.

- **Label elements**
- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms** GHS05, GHS07
- **Signal word** Danger
- **Hazard-determining components of labeling:**
guanidinium thiocyanate
- **Hazard statements**
Harmful if swallowed.
Causes severe skin burns and eye damage.
Harmful to aquatic life with long lasting effects.
- **Precautionary statements**
If medical advice is needed, have product container or label at hand.
Keep out of reach of children.
Read label before use.

(Contd. on page 2)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Binding Buffer

(Contd. of page 1)

- Do not breathe dusts or mists.
- Wash thoroughly after handling.
- Do not eat, drink or smoke when using this product.
- Avoid release to the environment.
- Wear protective gloves/protective clothing/eye protection/face protection.
- If swallowed: Call a poison center/doctor if you feel unwell.
- If swallowed: Rinse mouth. Do NOT induce vomiting.
- If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
- IF INHALED: Remove person to fresh air and keep comfortable for breathing.
- If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- Immediately call a poison center/doctor.
- Specific treatment (see on this label).
- Wash contaminated clothing before reuse.
- Store locked up.
- Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



Health = 3
Fire = 0
Reactivity = 0

· **HMIS-ratings (scale 0 - 4)**



Health = *3
Fire = 0
Reactivity = 0

· **Other hazards**

· **Results of PBT and vPvB assessment**

- **PBT:** Not applicable.
- **vPvB:** Not applicable.

3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** Mixture of the substances listed below with nonhazardous additions.

· **Dangerous components:**

CAS: 593-84-0	guanidinium thiocyanate	50%
CAS: 56-81-5	glycerol	50%

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Binding Buffer

(Contd. of page 2)

4 First-aid measures

· **Description of first aid measures**

· **General information:**

Immediately remove any clothing soiled by the product.

Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.

· **After inhalation:**

Supply fresh air. If required, provide artificial respiration if trained to do so. Keep patient warm. Consult doctor if symptoms persist.

In case of unconsciousness place patient stably in side position for transportation.

· **After skin contact:** Immediately wash with water and soap and rinse thoroughly.

· **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.

· **After swallowing:**

Rinse mouth

DO NOT induce vomiting.

· **Information for doctor:**

· **Most important symptoms and effects, both acute and delayed** No further relevant information available.

· **Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

5 Fire-fighting measures

· **Extinguishing media**

· **Suitable extinguishing agents:**

CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

Use fire fighting measures that suit the environment.

· **Special hazards arising from the substance or mixture**

Products of thermal decomposition of this material would include hydrogen cyanide, ammonia, and oxides of carbon, nitrogen and sulfur.

· **Advice for firefighters**

· **Protective equipment:** Wear self-contained breathing apparatus for fighting fires involving this material

6 Accidental release measures

· **Personal precautions, protective equipment and emergency procedures**

Wear self-contained breathing apparatus for responding to non-incident release of this material in which there is the potential for inhalation of vapors, mists or sprays

Wear protective equipment. Keep unprotected persons away.

· **Environmental precautions:**

Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

· **Methods and material for containment and cleaning up:**

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Use neutralizing agent.

Dispose contaminated material as waste according to item 13.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Binding Buffer

(Contd. of page 3)

Ensure adequate ventilation.

Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

Protective Action Criteria for Chemicals

PAC-1:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include guanidinium thiocyanate (0.98 mg/m³) and glycerol (45 mg/m³).

PAC-2:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include guanidinium thiocyanate (11 mg/m³) and glycerol (180 mg/m³).

PAC-3:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include guanidinium thiocyanate (65 mg/m³) and glycerol (1,100 mg/m³).

7 Handling and storage

Handling:

Precautions for safe handling

Ensure good ventilation/exhaustion at the workplace.

Prevent formation of aerosols.

Information about protection against explosions and fires: No special measures required.

Conditions for safe storage, including any incompatibilities

Store in cool, dry place. Store in well-ventilated location.

Storage:

Requirements to be met by storerooms and receptacles: No special requirements.

Information about storage in one common storage facility: Do not store together with acids or strong oxidizers

Further information about storage conditions: Keep receptacle tightly sealed.

Specific end use(s) Laboratory reagent

8 Exposure controls/personal protection

Additional information about design of technical systems: No further data; see item 7.

Control parameters

Work under a chemical fume hood when using this product. Ensure eyewash station and safety showers are readily accessible.

Components with limit values that require monitoring at the workplace:

The following constituent is the only constituent of the product which has a PEL, TLV or other recommended exposure limit.

At this time, the remaining constituent has no known exposure limits.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Binding Buffer

(Contd. of page 4)

CAS: 56-81-5 glycerol

Table with 2 columns: PEL, TLV. PEL: Long-term value: 15* 5** mg/m³ mist; *total dust **respirable fraction. TLV: TLV withdrawn-insufficient data human occup. exp.

Additional information: The lists that were valid during the creation were used as basis.

Exposure controls

The appropriate protective equipment under anticipated circumstances of use include lab-coat, safety glasses with side-shields and gloves.

Personal protective equipment:

General protective and hygienic measures:

- Keep away from foodstuffs, beverages and feed. Immediately remove all soiled and contaminated clothing. Wash hands before breaks and at the end of work. Avoid contact with the eyes and skin.

Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

Protection of hands:



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye protection:



Tightly sealed goggles



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Binding Buffer

(Contd. of page 5)

9 Physical and chemical properties

· Information on basic physical and chemical properties

· General Information

· Appearance:

Form:	Liquid
Color:	Clear
Odor:	Mild
Odor threshold:	Not determined.

· pH-value: Not determined.

· Change in condition

Melting point/Melting range:	Undetermined.
Boiling point/Boiling range:	Undetermined.

· Flash point: Not applicable.

· Flammability (solid, gaseous): Not applicable.

· Ignition temperature: 400 °C (752 °F)

· Decomposition temperature: Not determined.

· Auto igniting: Product is not selfigniting.

· Danger of explosion: Product does not present an explosion hazard.

· Explosion limits:

Lower:	0.9 Vol %
Upper:	0.0 Vol %

· Vapor pressure at 20 °C (68 °F): 0.1 hPa (0.1 mm Hg)

· Density: Not determined.

· Relative density Not determined.

· Vapor density Not determined.

· Evaporation rate Not determined.

· Solubility in / Miscibility with

Water: Fully miscible.

· Partition coefficient (n-octanol/water): Not determined.

· Viscosity:

Dynamic: Not determined.

Kinematic: Not determined.

· Solvent content:

Organic solvents: 50.0 %

VOC content: 0.00 %
0.0 g/l / 0.00 lb/gal

Solids content: 50.0 %

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Binding Buffer

(Contd. of page 6)

Other information No further relevant information available.

10 Stability and reactivity

- Reactivity No further relevant information available.
Chemical stability This product is normally stable under anticipated circumstances of use and storage.
Thermal decomposition / conditions to be avoided: Products of thermal decomposition of this material would include hydrogen cyanide, ammonia, and oxides of carbon nitrogen and sulfur.
Possibility of hazardous reactions No dangerous reactions known.
Conditions to avoid Avoid exposing product to extreme temperatures or incompatible chemicals
Incompatible materials: Acids and strong oxidizers
Hazardous decomposition products: Product will not undergo self-decomposition, so no such products will be generated.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:

LD/LC50 values that are relevant for classification:

CAS: 593-84-0 guanidinium thiocyanate

Oral LD50 593 mg/kg (rat)

- Primary irritant effect:
on the skin: Caustic effect on skin and mucous membranes.
on the eye: Strong caustic effect.
Sensitization: No sensitizing effects known.
Additional toxicological information:

The product shows the following dangers according to internally approved calculation methods for preparations:
Harmful
Corrosive
Swallowing will lead to a strong caustic effect on mouth and throat and to the danger of perforation of esophagus and stomach.

- Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

US

(Contd. on page 8)



Safety Data Sheet acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Binding Buffer

(Contd. of page 7)

12 Ecological information

Toxicity

Aquatic toxicity:

CAS: 593-84-0 guanidinium thiocyanate

EC50 42.4 mg/kg (daphnia)

Persistence and degradability No further relevant information available.

Behavior in environmental systems:

Bioaccumulative potential No further relevant information available.

Mobility in soil No further relevant information available.

Additional ecological information:

General notes:

Water hazard class 2 (Self-assessment): hazardous for water

Do not allow product to reach ground water, water course or sewage system.

Must not reach bodies of water or drainage ditch undiluted or unneutralized.

Danger to drinking water if even small quantities leak into the ground.

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

Other adverse effects No further relevant information available.

13 Disposal considerations

Waste treatment methods

Recommendation:

Dispose of contents in accordance with local/regional/national, and international recommendations.

Uncleaned packagings:

Recommendation:

Dispose of container in accordance with local/regional/national and international recommendations.

Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

UN-Number

DOT, IMDG, IATA

UN1760

UN proper shipping name

DOT

IMDG, IATA

Corrosive liquids, n.o.s. (guanidinium thiocyanate)

CORROSIVE LIQUID, N.O.S. (guanidinium thiocyanate)

(Contd. on page 9)



Safety Data Sheet acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Binding Buffer

(Contd. of page 8)

· Transport hazard class(es)

· DOT



· Class 8 Corrosive substances
· Label 8

· IMDG, IATA



· Class 8 Corrosive substances
· Label 8

· Packing group III
· DOT, IMDG, IATA

· Environmental hazards: Not applicable.

· Special precautions for user Warning: Corrosive substances
· Danger code (Kemler): 80
· EMS Number: F-A,S-B
· Stowage Category A
· Stowage Code SW2 Clear of living quarters.

· Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code Not applicable.

· Transport/Additional information:

· DOT
· Quantity limitations On passenger aircraft/rail: 5 L
On cargo aircraft only: 60 L

· IMDG
· Limited quantities (LQ) 5L
· Excepted quantities (EQ) Code: E1
Maximum net quantity per inner packaging: 30 ml
Maximum net quantity per outer packaging: 1000 ml

· UN "Model Regulation": UN 1760 CORROSIVE LIQUIDS, N.O.S. (GUANIDINIUM THIOCYANATE), 8, III

US

(Contd. on page 10)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Binding Buffer

(Contd. of page 9)

15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Sara**

· **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**

None of the ingredients is listed.

· **TSCA (Toxic Substances Control Act):**

All ingredients are listed.

· **Proposition 65**

· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms** GHS05, GHS07

· **Signal word** Danger

· **Hazard-determining components of labeling:**

guanidinium thiocyanate

· **Hazard statements**

Harmful if swallowed.

Causes severe skin burns and eye damage.

Harmful to aquatic life with long lasting effects.

· **Precautionary statements**

If medical advice is needed, have product container or label at hand.

Keep out of reach of children.

Read label before use.

Do not breathe dusts or mists.

Wash thoroughly after handling.

(Contd. on page 11)

**Safety Data Sheet**
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Binding Buffer

(Contd. of page 10)

Do not eat, drink or smoke when using this product.
Avoid release to the environment.
Wear protective gloves/protective clothing/eye protection/face protection.
If swallowed: Call a poison center/doctor if you feel unwell.
If swallowed: Rinse mouth. Do NOT induce vomiting.
If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
IF INHALED: Remove person to fresh air and keep comfortable for breathing.
If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
Continue rinsing.
Immediately call a poison center/doctor.
Specific treatment (see on this label).
Wash contaminated clothing before reuse.
Store locked up.
Dispose of contents/container in accordance with local/regional/national/international regulations.
• **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

• **Department issuing SDS:**

Zymo Research Corp.
Safety Department
17062 Murphy Ave.
Irvine, CA 92614
USA

Phone: 1-949-679-1190 or 1-888-882-9682

• **Contact:** sds@zymoresearch.com• **Date of preparation / last revision** 12/18/2018 / -• **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

Acute Tox. 4: Acute toxicity – Category 4

(Contd. on page 12)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 12/12

Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Binding Buffer

(Contd. of page 11)

Skin Corr. 1C: Skin corrosion/irritation – Category 1C

Eye Dam. 1: Serious eye damage/eye irritation – Category 1

Aquatic Chronic 3: Hazardous to the aquatic environment - long-term aquatic hazard – Category 3

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

1 Identification

- **Product identifier**
- **Trade name:** ZymoBIOMICS DNA Wash Buffer 1
- **Article number:** D4300-3-50, D4300-3-100
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**



GHS02 Flame

Flam. Liq. 3 H226 Flammable liquid and vapor.



GHS07

Eye Irrit. 2A H319 Causes serious eye irritation.

STOT SE 3 H336 May cause drowsiness or dizziness.

- **Label elements**

· **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms** GHS02, GHS07

· **Signal word** Warning

- **Hazard-determining components of labeling:**

propan-2-ol

ethanol

- **Hazard statements**

Flammable liquid and vapor.

Causes serious eye irritation.

May cause drowsiness or dizziness.

- **Precautionary statements**

If medical advice is needed, have product container or label at hand.

Keep out of reach of children.

Read label before use.

(Contd. on page 2)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 1

(Contd. of page 1)

Keep away from heat/sparks/open flames/hot surfaces. - No smoking.
Ground/bond container and receiving equipment.
Use explosion-proof electrical/ventilating/lighting/equipment.
Use only non-sparking tools.
Take precautionary measures against static discharge.
Avoid breathing dust/fume/gas/mist/vapors/spray
Wash thoroughly after handling.
Use only outdoors or in a well-ventilated area.
Wear protective gloves/protective clothing/eye protection/face protection.
If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
IF INHALED: Remove person to fresh air and keep comfortable for breathing.
If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
Continue rinsing.
Call a poison center/doctor if you feel unwell.
If eye irritation persists: Get medical advice/attention.
In case of fire: Use for extinction: CO2, powder or water spray.
Store in a well-ventilated place. Keep container tightly closed.
Store in a well-ventilated place. Keep cool.
Store locked up.
Dispose of contents/container in accordance with local/regional/national/international regulations.

Classification system:

NFPA ratings (scale 0 - 4)



HMIS-ratings (scale 0 - 4)



Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

3 Composition/information on ingredients

Chemical characterization: Mixtures

Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

Table with 2 columns: CAS number and percentage. Rows include propan-2-ol (25%) and ethanol (25%).

US

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 1

(Contd. of page 2)

4 First-aid measures

- **Description of first aid measures**
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:**
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:** Do not induce vomiting; immediately call for medical help.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed**
Inhalation of vapors, mists or sprays can cause drowsiness, dizziness, and other central nervous system effects. Accidental eye contact cause serious irritation.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures**
Wear protective equipment. Keep unprotected persons away.
- **Environmental precautions:** Dilute with plenty of water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Ensure adequate ventilation.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
- **Protective Action Criteria for Chemicals**

· **PAC-1:**

CAS: 67-63-0	propan-2-ol	400 ppm
CAS: 64-17-5	ethanol	1,800 ppm

· **PAC-2:**

CAS: 67-63-0	propan-2-ol	2000* ppm
CAS: 64-17-5	ethanol	3300* ppm

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 1

(Contd. of page 3)

· PAC-3:

CAS: 67-63-0	propan-2-ol	12000** ppm
CAS: 64-17-5	ethanol	15000* ppm

7 Handling and storage

- **Handling:**
- **Precautions for safe handling**
Ensure good ventilation/exhaustion at the workplace.
Prevent formation of aerosols.
- **Information about protection against explosions and fires:**
Keep ignition sources away - Do not smoke.
Protect against electrostatic charges.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep receptacle tightly sealed.
- **Specific end use(s)** Laboratory reagent

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**

· Components with limit values that require monitoring at the workplace:

CAS: 67-63-0 propan-2-ol	
PEL	Long-term value: 980 mg/m ³ , 400 ppm
REL	Short-term value: 1225 mg/m ³ , 500 ppm Long-term value: 980 mg/m ³ , 400 ppm
TLV	Short-term value: 984 mg/m ³ , 400 ppm Long-term value: 492 mg/m ³ , 200 ppm BEI
CAS: 64-17-5 ethanol	
PEL	Long-term value: 1900 mg/m ³ , 1000 ppm
REL	Long-term value: 1900 mg/m ³ , 1000 ppm
TLV	Short-term value: 1880 mg/m ³ , 1000 ppm

· Ingredients with biological limit values:

CAS: 67-63-0 propan-2-ol	
BEI	40 mg/L Medium: urine Time: end of shift at end of workweek Parameter: Acetone (background, nonspecific)

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 1

(Contd. of page 4)

· **Additional information:** The lists that were valid during the creation were used as basis.

· **Exposure controls**

· **Personal protective equipment:**

· **General protective and hygienic measures:**

- Keep away from foodstuffs, beverages and feed.
- Immediately remove all soiled and contaminated clothing.
- Wash hands before breaks and at the end of work.
- Avoid contact with the eyes.
- Avoid contact with the eyes and skin.

· **Breathing equipment:**

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

· **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:**



Tightly sealed goggles

9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· Form:	Liquid
· Color:	Clear
· Odor:	Alcohol-like
· Odor threshold:	Not determined.

· **pH-value:** Not determined.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 1

(Contd. of page 5)

· Change in condition	
Melting point/Melting range:	Undetermined.
Boiling point/Boiling range:	Undetermined.
· Flash point:	>25 °C (>77 °F)
· Flammability (solid, gaseous):	Not applicable.
· Ignition temperature:	425 °C (797 °F)
· Decomposition temperature:	Not determined.
· Auto igniting:	Product is not selfigniting.
· Danger of explosion:	Product is not explosive. However, formation of explosive air/vapor mixtures are possible.
· Explosion limits:	
Lower:	2 Vol %
Upper:	15 Vol %
· Vapor pressure at 20 °C (68 °F):	59 hPa (44.3 mm Hg)
· Density:	Not determined.
· Relative density	Not determined.
· Vapor density	Not determined.
· Evaporation rate	Not determined.
· Solubility in / Miscibility with Water:	Fully miscible.
· Partition coefficient (n-octanol/water):	Not determined.
· Viscosity:	
Dynamic:	Not determined.
Kinematic:	Not determined.
· Solvent content:	
Organic solvents:	50.0 %
VOC content:	50.00 %
	500.0 g/l / 4.17 lb/gal
Solids content:	0.0 %
· Other information	No further relevant information available.

10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:**
Products of thermal decomposition of this material would include carbon monoxide and carbon dioxide
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.

(Contd. on page 7)



Safety Data Sheet acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 1

(Contd. of page 6)

- Incompatible materials: No further relevant information available.
· Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects
· Acute toxicity:
· Primary irritant effect:
· on the skin: No irritant effect.
· on the eye: Irritating effect.
· Sensitization: No sensitizing effects known.
· Additional toxicological information:
The product shows the following dangers according to internally approved calculation methods for preparations:
Irritant

· Carcinogenic categories

· IARC (International Agency for Research on Cancer)

CAS: 67-63-0 propan-2-ol

3

· NTP (National Toxicology Program)

None of the ingredients is listed.

· OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
· Aquatic toxicity: No further relevant information available.
· Persistence and degradability No further relevant information available.
· Behavior in environmental systems:
· Bioaccumulative potential No further relevant information available.
· Mobility in soil No further relevant information available.
· Additional ecological information:
· General notes: Not hazardous for water.
· Results of PBT and vPvB assessment
· PBT: Not applicable.
· vPvB: Not applicable.
· Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
· Recommendation:
Dispose of contents in accordance with local/regional/national, and international recommendations.

(Contd. on page 8)



Safety Data Sheet acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 1

(Contd. of page 7)

- Uncleaned packagings:
Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Property and Value. Rows include UN-Number (UN1993), UN proper shipping name (Flammable liquids, n.o.s. (Ethanol, Isopropanol)), Transport hazard class(es) (3 Flammable liquids), Packing group (III), and Environmental hazards (Not applicable).

(Contd. on page 9)



Safety Data Sheet acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 1

(Contd. of page 8)

· Transport/Additional information:

· DOT

· Quantity limitations

On passenger aircraft/rail: 60 L
On cargo aircraft only: 220 L

· IMDG

· Limited quantities (LQ)

· Excepted quantities (EQ)

5L
Code: E1
Maximum net quantity per inner packaging: 30 ml
Maximum net quantity per outer packaging: 1000 ml

· UN "Model Regulation":

UN 1993 FLAMMABLE LIQUIDS, N.O.S. (ETHANOL, ISOPROPANOL), 3, III

15 Regulatory information

· Safety, health and environmental regulations/legislation specific for the substance or mixture

· Sara

· Section 355 (extremely hazardous substances):

None of the ingredients is listed.

· Section 313 (Specific toxic chemical listings):

CAS: 67-63-0 propan-2-ol

· TSCA (Toxic Substances Control Act):

All ingredients are listed.

· Proposition 65

· Chemicals known to cause cancer:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

· Carcinogenic categories

· EPA (Environmental Protection Agency)

None of the ingredients is listed.

· TLV (Threshold Limit Value established by ACGIH)

CAS: 67-63-0 propan-2-ol

A4

CAS: 64-17-5 ethanol

A3

(Contd. on page 10)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 1

(Contd. of page 9)

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

· **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms** GHS02, GHS07

· **Signal word** Warning

· **Hazard-determining components of labeling:**

propan-2-ol

ethanol

· **Hazard statements**

Flammable liquid and vapor.

Causes serious eye irritation.

May cause drowsiness or dizziness.

· **Precautionary statements**

If medical advice is needed, have product container or label at hand.

Keep out of reach of children.

Read label before use.

Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

Ground/bond container and receiving equipment.

Use explosion-proof electrical/ventilating/lighting/equipment.

Use only non-sparking tools.

Take precautionary measures against static discharge.

Avoid breathing dust/fume/gas/mist/vapors/spray

Wash thoroughly after handling.

Use only outdoors or in a well-ventilated area.

Wear protective gloves/protective clothing/eye protection/face protection.

If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

IF INHALED: Remove person to fresh air and keep comfortable for breathing.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

Call a poison center/doctor if you feel unwell.

If eye irritation persists: Get medical advice/attention.

In case of fire: Use for extinction: CO2, powder or water spray.

Store in a well-ventilated place. Keep container tightly closed.

Store in a well-ventilated place. Keep cool.

Store locked up.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

(Contd. on page 11)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 11/11

Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 1

(Contd. of page 10)

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 12/18/2018 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

BEI: Biological Exposure Limit

Flam. Liq. 3: Flammable liquids – Category 3

Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A

STOT SE 3: Specific target organ toxicity (single exposure) – Category 3

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

1 Identification

- **Product identifier**
- **Trade name:** ZymoBIOMICS DNA Wash Buffer 2
- **Article number:** D4300-4-60, D4300-4-200
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**



GHS02 Flame

Flam. Liq. 2 H225 Highly flammable liquid and vapor.



GHS07

Eye Irrit. 2A H319 Causes serious eye irritation.

STOT SE 3 H336 May cause drowsiness or dizziness.

- **Label elements**

- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).

- **Hazard pictograms** GHS02, GHS07

- **Signal word** Danger

- **Hazard-determining components of labeling:**

propan-2-ol

ethanol

- **Hazard statements**

Highly flammable liquid and vapor.

Causes serious eye irritation.

May cause drowsiness or dizziness.

- **Precautionary statements**

If medical advice is needed, have product container or label at hand.

Keep out of reach of children.

Read label before use.

(Contd. on page 2)



Safety Data Sheet acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 2

(Contd. of page 1)

Keep away from heat/sparks/open flames/hot surfaces. - No smoking. Ground/bond container and receiving equipment. Use explosion-proof electrical/ventilating/lighting/equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Avoid breathing dust/fume/gas/mist/vapors/spray Wash thoroughly after handling. Use only outdoors or in a well-ventilated area. Wear protective gloves/protective clothing/eye protection/face protection. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF INHALED: Remove person to fresh air and keep comfortable for breathing. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a poison center/doctor if you feel unwell. If eye irritation persists: Get medical advice/attention. In case of fire: Use for extinction: CO2, powder or water spray. Store in a well-ventilated place. Keep container tightly closed. Store in a well-ventilated place. Keep cool. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

Classification system:

NFPA ratings (scale 0 - 4)



HMIS-ratings (scale 0 - 4)



Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

3 Composition/information on ingredients

Chemical characterization: Mixtures

Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

Table with 2 columns: CAS number and percentage. Rows include CAS: 64-17-5 ethanol (25%) and CAS: 67-63-0 propan-2-ol (25%).



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 2

(Contd. of page 2)

4 First-aid measures

- **Description of first aid measures**
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:**
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:** Do not induce vomiting; immediately call for medical help.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures**
Wear protective equipment. Keep unprotected persons away.
- **Environmental precautions:**
Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Ensure adequate ventilation.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
- **Protective Action Criteria for Chemicals**

· **PAC-1:**

CAS: 64-17-5	ethanol	1,800 ppm
CAS: 67-63-0	propan-2-ol	400 ppm

· **PAC-2:**

CAS: 64-17-5	ethanol	3300* ppm
CAS: 67-63-0	propan-2-ol	2000* ppm

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 2

(Contd. of page 3)

PAC-3:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include ethanol (15000* ppm) and propan-2-ol (12000** ppm).

7 Handling and storage

- Handling:
Precautions for safe handling: Ensure good ventilation/exhaustion at the workplace.
Information about protection against explosions and fires: Keep ignition sources away - Do not smoke.
Conditions for safe storage, including any incompatibilities:
Storage:
Requirements to be met by storerooms and receptacles: Store in a cool location.
Information about storage in one common storage facility: Not required.
Further information about storage conditions: Keep receptacle tightly sealed.
Specific end use(s): Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
Control parameters

Components with limit values that require monitoring at the workplace:

Table with 2 columns: CAS number and exposure limits. Rows include ethanol (CAS: 64-17-5) and propan-2-ol (CAS: 67-63-0) with PEL, REL, TLV, and BEI values.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 2

(Contd. of page 4)

Ingredients with biological limit values:

CAS: 67-63-0 propan-2-ol

Table with 2 columns: BEI, and description. BEI: 40 mg/L. Medium: urine. Time: end of shift at end of workweek. Parameter: Acetone (background, nonspecific)

Additional information: The lists that were valid during the creation were used as basis.

Exposure controls

Personal protective equipment:

General protective and hygienic measures:

- Keep away from foodstuffs, beverages and feed. Immediately remove all soiled and contaminated clothing. Wash hands before breaks and at the end of work. Avoid contact with the eyes. Avoid contact with the eyes and skin.

Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

Protection of hands:



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye protection:



Tightly sealed goggles

US

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 2

(Contd. of page 5)

9 Physical and chemical properties

· Information on basic physical and chemical properties

· General Information

· Appearance:

Form:	Liquid
Color:	Clear
Odor:	Alcohol-like
Odor threshold:	Not determined.

· pH-value: Not determined.

· Change in condition

Melting point/Melting range:	Undetermined.
Boiling point/Boiling range:	78 °C (172.4 °F)

· Flash point: 13 °C (55.4 °F)

· Flammability (solid, gaseous): Not applicable.

· Ignition temperature: 425 °C (797 °F)

· Decomposition temperature: Not determined.

· Auto igniting: Product is not selfigniting.

· Danger of explosion: Product is not explosive. However, formation of explosive air/vapor mixtures are possible.

· Explosion limits:

Lower:	2 Vol %
Upper:	15 Vol %

· Vapor pressure at 20 °C (68 °F): 59 hPa (44.3 mm Hg)

· Density: Not determined.

· Relative density Not determined.

· Vapor density Not determined.

· Evaporation rate Not determined.

· Solubility in / Miscibility with

Water: Fully miscible.

· Partition coefficient (n-octanol/water): Not determined.

· Viscosity:

Dynamic: Not determined.

Kinematic: Not determined.

· Solvent content:

Organic solvents:	50.0 %
VOC content:	50.00 %
	500.0 g/l / 4.17 lb/gal

Solids content: 0.0 %

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 2

(Contd. of page 6)

Other information No further relevant information available.

10 Stability and reactivity

- Reactivity No further relevant information available.
Chemical stability
Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
Possibility of hazardous reactions No dangerous reactions known.
Conditions to avoid No further relevant information available.
Incompatible materials: No further relevant information available.
Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:
Primary irritant effect:
on the skin: No irritant effect.
on the eye: Irritating effect.
Sensitization: No sensitizing effects known.
Additional toxicological information:
The product shows the following dangers according to internally approved calculation methods for preparations:
Irritant

Carcinogenic categories

IARC (International Agency for Research on Cancer)

CAS: 67-63-0 propan-2-ol 3

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.
Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes:
Water hazard class 1 (Self-assessment): slightly hazardous for water

(Contd. on page 8)



Safety Data Sheet acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 2

(Contd. of page 7)

- Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation: Dispose of contents in accordance with local/regional/national, and international recommendations.
Uncleaned packagings:
Recommendation: Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Hazard/Requirement and Details. Rows include UN-Number (UN1993), UN proper shipping name (Flammable liquids, n.o.s.), Transport hazard class(es) (3 Flammable liquids), and IMDG, IATA details.

(Contd. on page 9)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 2

(Contd. of page 8)

Table with hazard information including Packing group (II), Environmental hazards (Not applicable), Special precautions (Warning: Flammable liquids), Transport information (Not applicable), and UN Model Regulation (UN 1993 FLAMMABLE LIQUIDS, N.O.S. (ISOPROPANOL, ETHANOL), 3, II).

15 Regulatory information

Table with regulatory information including safety regulations, Section 355 (None of the ingredients is listed), Section 313 (CAS: 67-63-0 propan-2-ol), TSCA (All ingredients are listed), and Proposition 65 (None of the ingredients is listed).

(Contd. on page 10)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 2

(Contd. of page 9)

Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

TLV (Threshold Limit Value established by ACGIH)

Table with 2 columns: CAS number and chemical name, and a third column for hazard codes (A3, A4). Rows include ethanol (CAS: 64-17-5) and propan-2-ol (CAS: 67-63-0).

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms GHS02, GHS07

Signal word Danger

Hazard-determining components of labeling:

propan-2-ol
ethanol

Hazard statements

Highly flammable liquid and vapor.
Causes serious eye irritation.
May cause drowsiness or dizziness.

Precautionary statements

If medical advice is needed, have product container or label at hand.
Keep out of reach of children.
Read label before use.
Keep away from heat/sparks/open flames/hot surfaces. - No smoking.
Ground/bond container and receiving equipment.
Use explosion-proof electrical/ventilating/lighting/equipment.
Use only non-sparking tools.
Take precautionary measures against static discharge.
Avoid breathing dust/fume/gas/mist/vapors/spray
Wash thoroughly after handling.
Use only outdoors or in a well-ventilated area.
Wear protective gloves/protective clothing/eye protection/face protection.
If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
IF INHALED: Remove person to fresh air and keep comfortable for breathing.
If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
Continue rinsing.
Call a poison center/doctor if you feel unwell.
If eye irritation persists: Get medical advice/attention.
In case of fire: Use for extinction: CO2, powder or water spray.
Store in a well-ventilated place. Keep container tightly closed.
Store in a well-ventilated place. Keep cool.
Store locked up.
Dispose of contents/container in accordance with local/regional/national/international regulations.

(Contd. on page 11)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 11/11

Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNA Wash Buffer 2

(Contd. of page 10)

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 12/18/2018 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

BEI: Biological Exposure Limit

Flam. Liq. 2: Flammable liquids – Category 2

Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A

STOT SE 3: Specific target organ toxicity (single exposure) – Category 3

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

1 Identification

- **Product identifier**
- **Trade name:** ZymoBIOMICS DNase/RNase Free Water
- **Article number:** D4302-5-1, D4302-5-3, D4302-5-10, D4032-5-50
- **CAS Number:**
7732-18-5
- **EC number:**
231-791-2
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The substance is not classified, according to the Globally Harmonized System (GHS).

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



Health = 0
Fire = 0
Reactivity = 0

- **HMIS-ratings (scale 0 - 4)**



Health = 0
Fire = 0
Reactivity = 0

- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.

(Contd. on page 2)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

· **vPvB:** Not applicable.

(Contd. of page 1)

3 Composition/information on ingredients

- **Chemical characterization: Substances**
- **CAS No. Description**
7732-18-5 water, distilled, conductivity or of similar purity
- **Identification number(s)**
- **EC number:** 231-791-2

4 First-aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.
- **After swallowing:** Do not induce vomiting; immediately call for medical help.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Use fire fighting measures that suit the environment.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Wear protective clothing.
- **Environmental precautions:** Dilute with plenty of water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 2)

· Protective Action Criteria for Chemicals

· PAC-1:

Substance is not listed.

· PAC-2:

Substance is not listed.

· PAC-3:

Substance is not listed.

7 Handling and storage

· Handling:

· Precautions for safe handling No special measures required.

· Information about protection against explosions and fires: No special measures required.

· Conditions for safe storage, including any incompatibilities

· Storage:

· Requirements to be met by storerooms and receptacles: No special requirements.

· Information about storage in one common storage facility: Not required.

· Further information about storage conditions: None.

· Specific end use(s) Laboratory reagent

8 Exposure controls/personal protection

· Additional information about design of technical systems: No further data; see item 7.

· Control parameters

· Components with limit values that require monitoring at the workplace: Not required.

· Additional information: The lists that were valid during the creation were used as basis.

· Exposure controls

· Personal protective equipment:

· General protective and hygienic measures:

The usual precautionary measures for handling chemicals should be followed.

· Breathing equipment: Not required.

· Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 3)

Eye protection: Goggles recommended during refilling.

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information

Appearance:

Form: Liquid
Color: Clear
Odor: Odorless
Odor threshold: Not determined.

pH-value: Not determined.

Change in condition

Melting point/Melting range: 0 °C (32 °F)
Boiling point/Boiling range: 100 °C (212 °F)

Flash point: Not applicable.

Flammability (solid, gaseous): Not applicable.

Decomposition temperature: Not determined.

Auto igniting: Not determined.

Danger of explosion: Product does not present an explosion hazard.

Explosion limits:

Lower: Not Applicable
Upper: Not Applicable

Vapor pressure at 20 °C (68 °F): 23 hPa (17.3 mm Hg)

Density at 20 °C (68 °F): 1 g/cm³ (8.345 lbs/gal)

Relative density: Not determined.

Vapor density: Not determined.

Evaporation rate: Not determined.

Solubility in / Miscibility with

Water: Fully miscible.

Partition coefficient (n-octanol/water): Not determined.

Viscosity:

Dynamic at 20 °C (68 °F): 0.952 mPas
Kinematic: Not determined.
Water: 100.0 %
VOC content: 0.00 %
0.0 g/l / 0.00 lb/gal

Solids content: 0.0 %

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 4)

Other information No further relevant information available.

10 Stability and reactivity

- Reactivity, Chemical stability, Thermal decomposition / conditions to be avoided, Possibility of hazardous reactions, Conditions to avoid, Incompatible materials, Hazardous decomposition products

11 Toxicological information

- Information on toxicological effects, Acute toxicity, Primary irritant effect, Sensitization, Additional toxicological information

IARC (International Agency for Research on Cancer)

Substance is not listed.

NTP (National Toxicology Program)

Substance is not listed.

OSHA-Ca (Occupational Safety & Health Administration)

Substance is not listed.

12 Ecological information

- Toxicity, Aquatic toxicity, Persistence and degradability, Behavior in environmental systems, Bioaccumulative potential, Mobility in soil, Additional ecological information, General notes

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 5)

- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:** Smaller quantities can be disposed of with household waste.
- **Uncleaned packagings:**
- **Recommendation:** Dispose of container in accordance with local/regional/national and international recommendations.
- **Recommended cleansing agent:** Water, if necessary with cleansing agents.

14 Transport information

· UN-Number	
· DOT, ADN, IMDG, IATA	not regulated
· UN proper shipping name	
· DOT, ADN, IMDG, IATA	not regulated
· Transport hazard class(es)	
· DOT, ADN, IMDG, IATA	
· Class	not regulated
· Packing group	
· DOT, IMDG, IATA	not regulated
· Environmental hazards:	Not applicable.
· Special precautions for user	Not applicable.
· Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code	Not applicable.
· UN "Model Regulation":	not regulated

15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Sara**

· Section 355 (extremely hazardous substances):
Substance is not listed.

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 6)

Section 313 (Specific toxic chemical listings):

Substance is not listed.

TSCA (Toxic Substances Control Act):

Substance is listed.

Proposition 65

Chemicals known to cause cancer:

Substance is not listed.

Chemicals known to cause reproductive toxicity for females:

Substance is not listed.

Chemicals known to cause reproductive toxicity for males:

Substance is not listed.

Chemicals known to cause developmental toxicity:

Substance is not listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

Substance is not listed.

TLV (Threshold Limit Value established by ACGIH)

Substance is not listed.

NIOSH-Ca (National Institute for Occupational Safety and Health)

Substance is not listed.

GHS label elements Void

Hazard pictograms Void

Signal word Void

Hazard statements Void

Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Department issuing SDS:

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

Contact: sds@zymoresearch.com

Date of preparation / last revision 12/18/2018 / -

Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

(Contd. on page 8)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 8/8

Safety Data Sheet
acc. to OSHA HCS

Printing date 12/18/2018

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 7)

IMDG: International Maritime Code for Dangerous Goods
DOT: US Department of Transportation
IATA: International Air Transport Association
ACGIH: American Conference of Governmental Industrial Hygienists
EINECS: European Inventory of Existing Commercial Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
NFPA: National Fire Protection Association (USA)
HMIS: Hazardous Materials Identification System (USA)
VOC: Volatile Organic Compounds (USA, EU)
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
NIOSH: National Institute for Occupational Safety
OSHA: Occupational Safety & Health
TLV: Threshold Limit Value
PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit

US