

5171 Wilfong Road Memphis, TN 38134

Telephone: 901-382-8716 Fax: 901-333-8223

Email: info@meridianlifescience.com www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note:	Centrifuge before	e opening to ensur	e complete recove	ery of vial contents.
1		· · · · · · · · · · · · · · · · · · ·	1	J

Catalog #: M01254M **Lot #:** 3G19115

Description: MAb to Human Alpha Fetoprotein

Monoclonal Antibody to Human Alpha Fetoprotein (AFP)

Specificity: Human AFP. Does not crossreact with Human Serum Albumin.

Clone: B491M

Isotype: IgG₁, kappa

Immunogen: Alpha Fetoprotein (AFP)

Format: Purified, Liquid

Purification: Ion Exchange Chromatography

Concentration: 6.66 mg/mL (OD280nm)

Affinity Constant: Not Determined

Buffer: PBS

Preservative: 0.09% Sodium Azide

Applications: Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its

particular application. Other applications have not been tested but use in such assays should not necessarily

be excluded.

Signature

Storage: Short term (up to 15 days) store at 2–8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

Warning: This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive

67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper

plumbing, flush with copious volumes of water to prevent azide build-up in drains.

25 Jan 2017

Date