



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: N86583H **Lot #:** 9A02916

Description: BNP and NT-proBNP Free Plasma
Pooled human plasma free of whole proBNP molecules and proBNP fragments (BNP and NT-proBNP).

Source: Pooled normal human plasma

Format: Affinity Purified, Liquid

Purification: Immunoaffinity chromatography. Gel for affinity column was prepared by coupling anti-BNP and anti-NT-proBNP monoclonal antibodies with BrCN-activated Sepharose CL 4B.

Concentration: Not Applicable

Buffer: Not Applicable

Preservative: None

Applications: Can be used for preparation of any proBNP-related standards and calibrators. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Store at -20°C. Avoid multiple freeze/thaw cycles.

Inactivation: Not Applicable

Warnings: All human source materials have tested negative for HIV1, HIV2, HCV antibodies and HBsAg. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.

Signature

29 Jan 2016

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY