



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134  
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)

[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** R01737 **Lot #:** 13K33217

**Description:** Enterovirus B Cox B3 Recomb.  
Enterovirus B Coxsackievirus B3 Recombinant  
Product contains a GST fusion partner.  
Molecular Weight: 60 kDa

**Source:** *E. coli*

**Format:** Purified, Liquid

**Purification:** > 85% Purity. Multi-step procedure including affinity chromatography.

**Concentration:** 1.04 mg/mL (OD280nm)

**Buffer:** Tris buffer with salts, pH 8.0.

**Preservative:** 0.09% Sodium Azide

**Applications:** Suitable for use as a control in Lateral Flow Assays and as an immunogen. Suitable for spiking into stool samples. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Storage:** Store at 2-8°C for < 1 week. For long term storage, aliquot and store at -20°C to avoid multiple freeze/thaw cycles.

**Inactivation:** Not Applicable

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

*Brenda Dunn*

29Nov2017

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**