



Meridian
Life Science,® Inc.
Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA
Telephone: 901-382-8716
Fax: 901-333-8223
Email: info@meridianlifescience.com
www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #:	R01742	Lot #:	4L33817
Description:	<i>E. coli</i> O157 VT2 Recombinant <i>Escherichia coli</i> O157 Verotoxin 2 (VT2) Recombinant Product contain the fusion partner His-tag, Molecular Weight: 36 kDa		
Source:	<i>E. coli</i>		
Format:	Purified, Liquid		
Purification:	Multi-step procedure including affinity chromatography.		
Concentration:	0.58 mg/mL (OD280nm)		
Buffer:	Carbonate Buffer with salts, pH 8.3.		
Preservative:	0.09% Sodium Azide		
Applications:	Suitable for use as a control in Lateral Flow Assays and as an immunogen. Suitable for spiking into stool samples. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage:	Store at 2-8°C.		
Inactivation:	Not Applicable		
Safety Note(s):	Refer to the appropriate Safety Data Sheet (SDS) for additional information.		

Diana Dunn

06Dec2017

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY