

CERTIFICATE OF ANALYSIS

Important Note:	Centrifuge before opening to ensure complete recovery of vial contents.			
Catalog #:	Y01273B		Lot #:	24L33817
Description:	THC BGG Conjugate Tetrahydrocannabinol (THC) Bovine Gamma Globulin Conjugate			
Source:	Not Applicable			
Format:	BGG, Liquid			
Concentration:	4.58 mg/mL			
Buffer:	Phosphate Buffered Saline, pH 7.4			
Preservative:	0.02% Sodium Azide			
Applications:	Suitable for use in ELISA and Lateral Flow assays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended antibody-antigen pair for ELISA and Lateral Flow assays: 			
Storage:	Short-term store at 2-8°C. Long term store at –20°C. Avoid multiple freeze/thaw cycles.			
Inactivation:	Not Applicable			
Safety Note(s):	Refer to the appropriate S	afety Data Sheet (SDS) for ad	ditional information.	

Gritcler Molsin

04 DEC 2017

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY